

A sentinel protein assay for simultaneously quantifying cellular processes



Martin Soste

Picotti Lab, ETH Zurich

Skyline user group meeting

ASMS 2015

perturbations →



ETH zürich



α -Synuclein induced cytotoxicity

1 copy α -Syn



2 copies α -Syn




➔ Established model for protein aggregation disease



collaboration with Lindquist lab

100 protective genes

Genetic modulator (1...100)



α -Syn + gene n





What are the genetic modulators doing to rescue the yeast?

How to find activated processes?

effects of 100 genes, triplicates = 300 samples


Cell biology & biochemistry



✓ functional activity

x time and labour


Omics



✓ unbiased

x not necessarily activity

Targeted proteomics



✓ quantitative

x biased

x not necessarily activity


Methodological aim

- unbiased screening for activity with targeted proteomics

✓ unbiased


✓ functional activity

✓ quantitative




The sentinel protein assay

- A fingerprint for activity




Effectors of protective genes?




Sentinels:

Markers that report cellular process activation


Sentinels are proteins, phosphorylation sites or degradation products



Example: phosphorylation-based sentinels




Example: MAPK activation



Identified functionally validated sentinels

SENTINEL
ATG1			
ATG13			
IRE1			
HOG1			
FUS3			
SLT2			
CK1			
MAF1			
SCH9			
YPK1			
IGO1			
MIG1			
SUI2			
DUN1			
SLX4			
TLG1			
TLG1			
SIP1			
CTK1			
ADR1			
SEC63			
CDC23			
CDC24			
CDH1			
CDH1	P-site	CDC28 - KINASE ACTIVITY	Phosphorylated by CDC28, in S, G2 and M phase
Bcy1	P-site	YAK1 - KINASE ACTIVITY	Phosphorylated by YAK1 in response to glucose star
YPK1	P-site	PKH1 - KINASE ACTIVITY	Phosphorylated by Pkh1
NPR1	P-site	NITROGEN METABOLISM	Dephosphorylated upon nitrogen limitation
RCK2	P-site	HOG1 - KINASE ACTIVITY	Phosphorylated by HOG1 after osmotic stress.
(.....)		(.....)	(.....)



Fingerprint assay:

- 157 proteins
- 152 P-sites

SRM assays



- Proteotypic peptides
PROT – 300, PHOS – 152
- Top 3–6 transitions
PROT – 1174, PHOS – 807
- iRT-based scheduling
- Synthetic heavy peptide confirmation



188 processes

Testing the sentinel protein assay

Well-characterized perturbations

Heat shock

Heat-shock recov.


Rapamycin

AA/N starvation


Osmo-stress

Osmo-stress adapt.

Quiescence









“Sentinel Fingerprints”




Multi-peptide chromatograms

Validation from expected responses




Simultaneous and large-scale coverage




- Snapshots of yeast physiology
- Report of 202 markers for activity
- Success rate = 74%

Novel responses in stationary phase




Concentration of cytosol



- Water deficiency marker up

Coordinated response of HSPs




- Subset of heat-shock proteins up
- Typical markers of heat-shock down

Summary: sentinel assay

Advantages:

- Captures behaviour of many pathways
- Probes activation, not merely abundance
- Targeted yet unbiased (“discovery SRM”)
- Fast and interpretable readout of cell status




➔ Apply to any perturbation of yeast
(e.g. large set of genetic manipulations or drugs)

Limits:


- Cannot identify a perturbed pathway if a marker is not included
- Influence of pathway crosstalk

Outlook: fingerprints of rescue from α -Syn



Protective genes (~100)



Sentinels



Multivariate regression



- Identify key processes in rescue
- Modulate key processes in higher organism models of α -Syn toxicity

Acknowledgements

Supervisor

Paola Picotti

Picotti Lab

Andre Melnik

Yuehan Feng

Juan Gerez

Rita Hrabakova

Paul Boersema


Pascal Leuenberger

Ilaria Piazza

Simone Schopper

Natalia Prymaczok

Marco Tognetti



α -Syn model

Susan Lindquist, WIBR

Gabriela Caraveo, WIBR

Sentinel Prediction

Stefanie Wanka, UZH

Christian Von Mering, UZH

Data analysis

Andreas Beyer, Sybacol

Picture on title slide
from Cooper *et al.* 2006



Thank you to Brendan, the development team and the user community!

ETH zürich



FNSNF



European
Research
Council



martin.soste@bc.biol.ethz.ch