



Skyline Tutorial Webinar #6

Welcome to the sixth in a series of tutorial webinars designed to help you get the most out of Skyline proteomics software.

We'll begin shortly ...

▶ **Webinar Audio:** Make sure you have audio during the webinar.

Click on the telephone icon in the JoinMe toolbar - **circled** here →

You then have the choice to either:

- ▶ Connect by phone (US and International numbers)
- ▶ Connect by internet – a short, temporary software download is required.



▶ **Submit Skyline Questions:** Use this form to submit Skyline questions to the team – we will try to answer as many as possible:

- ▶ The QA link will also be available in the chat window – to open your chat window click the bubble icon – **circled** here →





Skyline Tutorial Webinar #6

Effective Data Processing and Interrogation with Skyline


With
Brendan MacLean (Principal Developer)

Agenda


- ▶ Welcome from the Skyline team!
- ▶ **Effective Data Processing and Interrogation**
- ▶ Introduction with Brendan MacLean
 - ▶ Brief overview of refinement by peptide/protein differentiation
- ▶ Tutorial with Brendan MacLean
 - ▶ Effective grouped study data processing and QC
 - ▶ Sneak peak: Group Comparison features in Skyline 3.1
- ▶ Audience Q&A – submit questions to Google Form:
<http://tinyurl.com/QA4Skyline>



Supporting a Broader Hypothesis



Targeted Method Refinement



Identifying Plasma Proteins with Altered Levels in a Dahl Salt Sensitive Rat Model

7 Rats on a Low Salt Diet




7 Rats on a High Salt Diet



109 proteins reported in literature as related to heart disease




Phases of Refinement: Detection



Scheduling Multi-Replicate Assays

- ▶ 49 Proteins
- ▶ 135 Peptides
- ▶ 780 Transitions
- ▶ 6 Unscheduled Injections



Prior Knowledge and Consistency



- ▶ Based on empirical measurement
- ▶ Powerful enough to be used cross-lab / cross experiment
- ▶ More powerful run-to-run

- ▶ Relative ion abundance
 - ▶ Spectral and chromatogram libraries
- ▶ Retention time
 - ▶ iRT
- ▶ Optimal collision energy
- ▶ Ion mobility
 - ▶ Collisional cross-section




Haptoglobin

LQTEGDGIYTLNSEK




Haptoglobin

LQTEGDGIYTLNSEK





SVVDIGLIK




Mitochondrial 39S ribosomal protein L9

CSSLWAGAAWLR




Mitochondrial 39S ribosomal protein L9

CSSLLWAGAAWLR




SVVDIGLIK







CSSLWAGAAWLR

H_159_REP1





D_138_REP3



Truncated and Missing Peaks

TGTLNLMDFLSR



Tutorial

- ▶ Effective Data Processing and Interrogation with Skyline



Learn More

- ▶ **NEW!** Grouped Studies Tutorial
- ▶ Webinar #7: iRT or Isotope Labeled Standards
 - ▶ Tuesday, May 12
- ▶ Workshop and ASMS
- ▶ Skyline User Group Meeting at ASMS
 - ▶ May 31 at Old Post Office, St. Louis, MO
- ▶ Workshop in Rio de Janeiro, August 31-September 2
- ▶ Workshop in Puerto Vallarta, November
- ▶ Weeklong Course at PRBB, Barcelona,
 - ▶ November 15-20



Questions?

- ▶ Ask any questions you have on this month's topic, Effective Data Processing, at the following form:

<http://tinyurl.com/QA4Skyline>

- ▶ Take the post-webinar survey

<http://tinyurl.com/Survey4Webinar>





Skyline Tutorial Webinar #6

This ends this Skyline Tutorial Webinar.

Please give us feedback on the webinar at the following survey:

<http://tinyurl.com/Survey4Webinar>

A recording of today's meeting will be available shortly at the Skyline website.

We look forward to seeing you at a future Skyline Tutorial Webinar.

